

PRACTICE TIPS

How to use the barcode labels correctly

Only
1 barcode set
per patient

IMPORTANT

Use a new barcode set for each patient.

- The labels can be used in any order.
- You can order the barcode labels using our order form.

This barcode is intended for the **request form**.

Schein


10 1000 1258 01	10 1000 1258 02	10 1000 1258 03
10 1000 1258 04	10 1000 1258 05	10 1000 1258 06
10 1000 1258 07	10 1000 1258 08	10 1000 1258 09
10 1000 1258 10	10 1000 1258 11	10 1000 1258 12
10 1000 1258 13	10 1000 1258 14	10 1000 1258 15

These labels are provided for the **patient's sample material**.

These **number codes** are intended for small sample tubes. (e.g. CapiSave)

This label is provided for your **laboratory journal**.

10 1000 1258	10 1000 1258	10 1000 1258
Name:		
Date of birth:		
10 1000 1258		



DFÜ order number

xx XXXX XXXX 00

10 1000 1258 01

10 1000 1258 01

You can enter this DFÜ order number in your IT system or scan the barcode.



biovis Diagnostik MVZ GmbH
Brüsseler Str. 18 | 65552 Limburg-Eschhofen
info@biovis.de | biovis.de

